

Introduction

- Capsule filling processes via a dosator are widely applied in the pharmaceutical industry. Capsule filling of **spray dried powders** using a dosator-based capsule filler can be challenging due to the **cohesive properties** inherent to those powders.

MG2[®] Flexalab is a dosator-based capsule filler suitable for precision capsule filling, integrated with a 100% weight control system, MultiNETT, controlling in process the net weight contained in each single capsule.

- For carrier-based powders, two main attributes were identified as major players in a low-dosage dosator-based capsule filling process: **the ratio between the dosing chamber length and powder layer height** and a **homogenous powder layer**^[1, 2].

The **main goal** of this work were to assess precision capsule filling of a model spray dried powder using a dosator-based **MG2[®] Flexalab unit**, optimize the filling process and evaluate its impact on powder in-vitro aerodynamic performance.

Materials and Methods

Capsule filling set-up



Figure 1 – MG2[®] Flexalab set-up: machine, rotary container and dosator.

Dosator-based filling mechanism




Figure 2 – Schematic representation of a dosator-based filling mechanism.

Spray drying process parameters

Table 1 – Spray drying process parameters.

Composition (% w/w)	Solids concentration (%)	Solvent system (% w/w)	Feed rate (g/min)	Atomizing rate (mm in rotameter)	Dry gas flow rate (kg/h)	Outlet temperature (°C)
Trehalose: L-leucine 80:20	2	Water: Ethanol 50:50	7	50	35	70

Experimental Design

Capsule	Experimental Set-up	Aerodynamic Performance
Clear HPMC Size #3	Speed (caps/h) 1500 2000	Andersen Cascade Impactor
Fill weigh	Dosator Diameters (mm) 1.9 2.2 2.8 3.4 3.7	Plastiapipe 60 L/min (4 kPa pressure drop)
20 mg ±1 10 mg ±1 5 mg ±0.5	Chamber / Layer ratio From 0.5 To 0.9	N=3

Results and Discussion

Capsule filling optimization

- Powder adhesion:** internal walls and behind the scrapper
- Agglomerates** in powder layer
- Powder accumulation** around the dosator
- High powder compaction** in layer
- High rejection rates** during capsule filling



Figure 3 – Set-up and powder appearance before optimization.

Powder characterization




Figure 4 – SEM images of the spray dried powder.

- Cleaning system** for the removal of the excess powder from the dosator;
- Engineered mixing powder rod** to homogenize the powder bed and decrease the powder adhesion to the walls of the rotary container;

Table 2 – SD powders capsule filling optimization: process parameters and results for five experimental runs.

Run	Fill weight (mg)	Rejection rate (%)	Dosator diameter (mm)	Chamber / Layer ratio	Speed (caps./h)	Capsule visual observation
1	20	0	2.8	0.9	1500	High
2	20	0	3.4	0.7	1500	Medium
3	20	3	3.7	0.5	2000	Low
4	10	0	2.8	0.6	2000	Low
5	5	0	1.9	0.8	2000	Low

Aerodynamic performance




Figure 5 – Aerodynamic performance results measured by ACI for the different fill weights.

Aerodynamic performance in accordance with **visual observation** of capsules.

Multivariate statistical analysis




Figure 6 – Regression coefficient plots for ED_{Lc} (%) and FPD (µg/capsule)

Conclusion

Capsule filling process of spray dried composite particles using a MG2 Flexalab machine was **successfully achieved**.

Low powder compaction in capsules and low rejection rates were possible to obtain by **optimizing process parameters** and by **implementing appropriate engineering solutions**.

Good aerodynamic performances were obtained using a **reliable and robust technology** in a manufacturing environment, which is **easily scaled-up**.